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INTENDED USE	SAMPLES TESTED	POSITIVE RESULTS	PERFORMANCE	STUDY
Trisomy 21	1,696	210 of 212	Sensitivity: 99.1% Specificity: 99.9%	DNA sequencing of maternal plasma to detect Down syndrome: an international clinical validation study. Genet Med. 2011;13(11):913-920.
Trisomy 18	1,988	59 of 59	Sensitivity: >99.9% Specificity: 99.6%	DNA sequencing of maternal plasma identifies trisomy 18 and trisomy 13, as well as Down syndrome: an international collaborative study Genet Med. 2012;14(3):296-305.
Trisomy 13		11 of 12	Sensitivity: 91.7% Specificity: 99.7%	
Multiple gestations	25	8 of 8 trisomies (7 trisomy 21, 1 trisomy 13)	Detection rate: >99.9%	DNA sequencing of maternal plasma to identify Down syndrome and other trisomies in multiple gestations. Prenat Diagn. 2012;32(8):730-734.
INTENDED USE	SAMPLES STUDIED		PERFORMANCE	STUDY
Y chromosome	2,107		Accuracy: 99.4%	Accuracy of Noninvasive Prenatal Sex Determination using Massively Parallel Sequencing in Samples from a Large Clinical Validation Study. Mazloom AR, et al. American Society of Human Genetics. November 2012.
COMBINED SEX ANEUPLOIDIES	SAMPLES TESTED	POSITIVE RESULTS	PERFORMANCE	SOURCE
45,X; 47,XXX; 47,XXY; 47,XYY	420	25 of 26	Sensitivity: 96.2% Specificity: 99.7%	Noninvasive prenatal detection of sex chromosomal aneuploidies by sequencing circulating cell-free DNA from maternal plasma. Mazloom AR, et al. Prenat Diagn. 2013 Jun;33(6):591-597.

Clinical outcome of subchromosomal events detected by whole-genome noninvasive prenatal testing.

Helgeson J et al., Prenat Diagn. 2015 Oct;35(10):999-1004. doi: 10.1002/pd.4640. Epub 2015 Jul 27.

Table 1 Clinical outcomes of noninvasive prenatal testing identified deletions

Deletion; associated syndrome (n)	Total identified	Confirmed true positive	Suspected due to clinical findings	No additional information	Confirmed false positive	PPV lower-upper estimate (95% confidence) (%)
22q11.2; DiGeorge (175 393)	32	23	8	1	0	96.9-100 (82.0-100)
1p36- (175 393)	5	3	0	1	1	60.0-80.0 (17.0-98.9)
15q; Prader-Willi/Angelman (175 393)	9	8	1	0	0	100.0 (59.8-100.0)
5p; Cri-du-chat (175 393)	6	4	0	0	2	66.7 (24.1-94.0)
4p; Wolf-Hirschhorn (52 297)	1	1	0	0	0	100.0 (5.5-100.0)
11q; Jacobsen (52 297)	1	1	0	0	0	100.0 (5.5-100.0)
8q; Langer-Giedion (52 297)	1	0	1	0	0	100.0 (5.5 - 100.0)

Confirmed cases were determined by diagnostic testing (chromosomal microarray analysis, karyotyping, or fluorescence *in situ* hybridization) of specimens obtained by an invasive procedure in fetal, maternal, or both. Suspected cases were those where diagnostic testing was declined and phenotypic data via clinical presentation in person or by ultrasound were consistent with the deletion. Positive predictive value lower estimate is calculated presuming all patients with no additional information are confirmed falsely positive. Positive predictive value (PPV) upper estimate is calculated presuming all patients with no additional information are confirmed true positives. The 95% confidence level is calculated according to the efficient-score method and corrected for continuity.